



Protocol Calretinin - ELISA

Products

30% H₂O₂ solution: Merck, 8222887.1000
3,3',5,5'-tetramethylbenzidin (TMB): Fluka, 87748
Aceton: Fluka, 00570
Bovine serum albumin: Fluka, 05480
Citric acid monohydrate: Fluka, 27490
Ethanol: Fluka, 02860
Kathon MW/WT: Christ AG, Aesch, Switzerland
NaHCO₃: Fluka, 71628
Phenol: Merck, 206
Sulfuric acid 97%: Fluka, 84720
Tris/HCl (Tris-(hydroxymethyl)-aminoethan): Fluka, 93352
Tween-20: Fluka, 93773

Reagents

REAGENT 1. Rabbit CR 7679 antibody solution

- Tris/HCl pH 7.5 20 mmol/L
- NaCl 150 mmol/L
- Kathon MW/WT 0.2 ml/L
- Amaranth 10 mg/L
- **Added only before use:** rabbit CR 7697 antibody (1:5)

This antibody solution is used for sensitization of micro-well plate. Preparation: CR 7697 antibody is added ONLY directly before running the test at a concentration of 1:5! REAGENT 1 will be diluted at 1:1'000 in SOLUTION #3.1. -> end concentration of the antibody 1:5'000.

Stability: Kathon allows this solution to be stored at least 6 month at 2-8°C without significant loss of immunological activity

REAGENT 2. Mouse CR 6B3 antibody solution

- Tris/HCl pH 7.5 0.1 mol/L
- Bovine serum albumin (BSA) 5.0 g/L
- Phenol 1.0 g/L
- Kathon (MW/WT) 0.2 ml/L
- Evans Blue 10 mg/L
- **Added only before use:** mouse CR 6B3 antibody (1:10 - 1:20)

Preparation: monoclonal mouse CR 6B3 antibody is added ONLY directly before running the test at a concentration of 1:10! REAGENT 2 will be diluted 1:500 in REAGENT 4 -> end concentration of the antibody 1:5'000.

Stability: at least 6 months at 2-8°C



REAGENT 3. Calretinin standard solution

1000 ng/mL recombinant rat Calretinin in:

- Potassium-phosphate pH 6.5 0.1 mol/L
- Bovine serum albumin (BSA) 5.0 g/L
- Kathon MW/WT 0.2 ml/L

The recombinant Calretinin is ONLY added to the solution directly before use at a concentration of 1ug/ml!

Stability: at least 6 months at 2-8°C

REAGENT 4. Test dilution buffer

- Tris/HCl pH 7.5 0.2 mol/L
- Bovine serum albumin (BSA) 5.0 g/L
- Kathon MW/WT 0.2 ml/L
- Phenol 1.0 g/L

Preparation: add 6.05 g of Tris/HCl, 0.1 ml kathon MW/WT, 0.5 g phenol and 2.5 BSA in 400 ml ddH₂O. Adjust pH to 7.5 with 4 N HCl, and complete the volume to 500 ml with ddH₂O.

Stability: at least 6 months at 2-8°C

REAGENT 5. Tetramethylbenzidin-H₂O₂ solution (TMB)

- 3,3',5,5'-tetramethylbenzidin 20 mmol/L
- Aceton 100 ml/L
- Ethanol 900 ml/L
- H₂O₂ 50 mmol/L

Preparation: dissolve 240 mg TMB in 5 ml aceton. Further add 45 ml ethanol and then 300 µl of 30% H₂O₂ solution.

Stability: at least 6 month at 15-25°, in a well closed, light-protected glass container

REAGENT 6. Substrate buffer

- Potassium-citrate buffer pH 4.1 30 mmol/L
- Kathon MW/WT 0.2 ml/L

Preparation: dissolve 6.3 g citric acid monohydrate in about 800 ml ddH₂O.

Adjust pH to 4.1 with 4 N KOH. Complete the volume to 1000 ml with ddH₂O and add 200 µl kathon MW/WT

Stability: Few months at 15-25°C

Other solutions

Sensitization buffer (SOLUTION #3.1.)

- NaHCO₃ pH 9.5 100 mmol/L

Preparation: dissolve 4.2 g NaHCO₃ in 400 ml ddH₂O. Adjust pH to 9.5 with 4 N NaOH. Complete the volume to 500 ml with ddH₂O

Stability: at least 6 months at 2-8°C



Blocking buffer (SOLUTION #3.2.)

- Tris/HCl pH 7.5 200 mmol/L
- Bovine serum albumin (BSA) 10 g/L
- Kathon MW/WT 0.2 ml/L

Preparation: dissolve 12.1 g of Tris-HCl and 0.1 ml Kathon MW/WT in 400 ml ddH₂O. Adjust pH to 7.5 with 4 N HCl. Further add 5 g BSA. Complete the volume to 500 ml with ddH₂O

Stability: at least 6 months at 2-8°C

Stop buffer (SOLUTION #3.3.)

- Sulfuric acid 1.0 mol/L

Preparation: add 28 ml sulfuric acid to 500 ml ddH₂O

→ Warning: do not add water in acid

Washing buffer (SOLUTION #3.4.)

- Tween-20 0.5 ml/L

Preparation: dilute 500 µl Tween-20 in 1000 ml ddH₂O

Other materials

Micro-well plates: polystyrol, flat bottom, 96-wells, NUNC-immunoplate, MAXISORP, 4-42 404

Multichannel photometer: filter 450 nm

ELISA running

Sensitization of micro-well plate with rabbit CR 7679 antibody

20 ml sensitization solution is required per micro-well plate

- A corresponding amount of antibody solution (**REAGENT 1**) is diluted 1:1000 in sensitization buffer (**SOLUTION #3.1.**): for example 20 µl **REAGENT 1** in 20 ml **SOLUTION #3.1.**
→ rabbit CR 7679 antibody (end concentration 1:5'000) in 100 mmol/L NaHCO₃ pH 9.5
- Pipette 200 µl of this solution per well
- Incubate 16-24 hrs at 15-25°C
- Remove sensitization solution and wash 3x with ddH₂O (300 µl/well)
- Pipette 200 µl of blocking buffer (**SOLUTION #3.2.**) per well
- Close the plate with parafilm and incubate at least one day at 15-25°C

The micro-well plates sensitized with CR 7679 antibody can be stored (under humidified conditions at 2-8°C) in the blocking buffer until 2 months without significant loss of immunological binding capacity

Preparation of solutions for standard curve and loading of samples

- Dilute Calretinin standard solution 1000 ng/ml (**REAGENT 3**) in test dilution buffer (**REAGENT 4**)
→ 1:200 for 5.0 ng/ml
→ 1:400 for 2.5 ng/L



→ 1:800 for 1.25 ng/L

→ The null-value is given by the test dilution buffer alone

- h) Prepare mouse CR 6B3 antibody test solution: Dilute mouse CR 6B3 antibody solution (**REAGENT 2**) 1:500 in test dilution buffer (**REAGENT 4**) (mouse CR 6B3 antibody end concentration 1:5'000)

5 ml of this preparation is required per micro-well plate (5 ml **REAGENT 4** and 10 µl **REAGENT 2**)

- i) Remove blocking buffer from the wells
j) Pipette 200 µl of standard curve solutions and samples (user-defined dilution with **REAGENT 4**) in duplicates in each well, and add 50 µl of CR 6B3 antibody test solution prepared in step h)
k) Close the plate and incubate it for 16-24 hrs at 2-8°C (possible to do it over the week end)
l) Empty the wells and wash 3x with washing buffer (**SOLUTION #3.4.**) (300µl/well)
m) Incubation with 200 µl of rabbit anti-mouse-HRP (DAKO P0260, 1:1000 in **REAGENT 4**) for 4-5 hrs at room temperature.
n) Empty the wells and wash 3x with washing buffer (**SOLUTION #3.4.**) (300µl/well)
o) Wash 2-3x with dd. H₂O (300µl/well).

Enzyme reaction

- p) Mix 1 volume of TMB-H₂O₂ solution (**REAGENT 5**) with 20 volumes of substrate buffer (**REAGENT 6**)

21 ml of this solution is required per micro-well plate

This solution should not be vortexed (could cause precipitation) and has to be prepared just before use (stability 1 hr at 15-25°C)

- q) Pipette 200 µl of this TMB-H₂O₂ solution per well
r) Incubate 10 minutes at 15-25°C
s) Stop the reaction by adding 100 µl stop buffer (**SOLUTION #3.3.**)
t) Within 1 hr, analyze the plate with the photometer at 450 nm